
Cadwell® Sleepmate

Instructions for Use (IFU)

Pressure Transducer

| Contents | Page |
|-----------------|-------------|
| English | 4 |
| Dansk..... | 6 |
| Deutsch..... | 8 |
| Español | 10 |
| Suomi | 12 |
| Français..... | 14 |
| Italiano..... | 16 |
| Nederlands..... | 18 |
| Norsk | 20 |
| Português..... | 22 |
| Svenska..... | 24 |

Intended use

The Cadwell Sleepmate™ AC pressure airflow sensor is to detect respiratory airflow and/or snoring for recording onto a polygraph, EEG or similar type of physiological recorder. The sensor plugs directly into the electrode jack box; no external power source is required. A nasal cannula is attached to the patient.

Instructions for use

Attach either an ET_{CO}₂ or standard O₂ cannula to the sensor and to the patient to detect the airflow. The airflow sensed through the nose will then provide both a nasal airflow signal and snore signal. There are two sets of cables for connection to the recorder. One is for the airflow channel and one is for the snore channel. The two black connectors are connected internally in the sensor. The sensor can be connected to a nasal CPAP mask using standard respiratory tubing instead of using a cannula. When connecting an O₂ cannula to the sensor always use the Luer® connector to connect to the sensor box. The ET_{CO}₂ cannula should be connected directly to the sensor box. The recorder's High Pass filter should be set at the lowest possible setting in order to maximize the time constant.

Ensure the sensor is securely attached to the patient as described and you have clear, reliable signals prior to lights out.

Cleaning

This sensor is reusable and can be cleaned by rubbing with any cleansing agent that is **not** corrosive to plastic. The cannula used is disposable and should be changed between patients. For sterilization implement a standard procedure for gas sterilization. Do not autoclave this product. Use the molded ends of the jackpins when removing the connector from the junction box.

Technical Specifications

| Description | Properties |
|-------------------|--|
| Lead wire | Separate output cables for flow and snore channels |
| Connector | 1.5 mm (0.1") safety DIN connectors |
| Sensor | Piezo |
| Cannula Connector | Luer® |
| Filter Settings | Flow: High Pass: 0.01 Hz; Low Pass: 30 Hz Snore: High Pass: 30 Hz; Low Pass: NA |

Caution Statement:

US Federal Law restricts this device to sale by or on the order of a physician. This product is for diagnostic purposes only and is not intended to be used as an apnea monitor or in a life supporting or sustaining capacity.

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©2021 Cadwell Industries, Inc. All rights reserved. PN 302993-639-1 Rev. 00
Cadwell Industries Inc.
1 (800) 245-3001 U.S.A. toll-free | +1 (509) 735-6481
This document contains trademarks that belong to Cadwell Industries, Inc.