
Cadwell® Sleepmate

Instructions for Use (IFU)

CannuTherm™ Airflow Sensor

Contents	Page
English	4
Dansk.....	6
Deutsch.....	8
Español.....	11
Suomi	14
Français.....	16
Italiano.....	19
Nederlands.....	22
Norsk.....	24
Português.....	26
Svenska.....	28

Intended use

The Cadwell Sleepmate™ Cannutherm™ thermal/pressure airflow sensor is intended to detect patient nasal and/or oral breathing during sleep research and clinical studies of sleep.

Instructions for use

The thermal sensor should be inserted into the port on the cannula so that the short nasal probe is inside the cannula chamber where it is exposed to the nasal airflow. The oral sensor probe is flexible to allow it to be bent **gently** to insure that it is in front of the patient's mouth and in the path of oral airflow. The oral probe should not touch the skin as contact may attenuate the signal amplitude. Once the thermal sensor is inserted in the CannuTherm™ cannula, the lead wire should be clipped into the clips provided on the cannula tube. The entire sensor should be placed onto the patient in a manner similar to that of an oxygen cannula with the thermal lead wire and the cannula tubes looping over the ears.

Connect the sensor to the appropriate channel input of your recorder.

Ensure the sensor is securely attached to the patient as described and you have clear, reliable signals prior to lights out.

The sensors are delicate and must be handled carefully. Proper handling and cleaning will increase durability and longevity in the life of the sensor. The assembly should be carefully secured using both the sliding cinch tubing and surgical tape.

Cleaning

The thermal sensor was designed to be reusable and can be cleaned by wiping with any cleansing agent that is **not** corrosive toward plastic. **The entire assembly must be completely dry before using again to insure patient safety.** For sterilization implement a standard procedure for gas sterilization. Do not autoclave this product. Use the molded ends of the jackpins when removing the connector from the junction box. The cannula portion is single use and should be disposed of properly after use.

Technical Specifications

Description	Properties
Lead wire	Two conductor white lead wire, approximately 200 cm long.
Sensor	Flexible, coated probe for oral airflow, epoxy coated nasal airflow probe.
Connectors	Two Safety pins 1.5 mm (0.1"). Some recorders may use Lemo or other connectors.
Time constant	One second or longer. The use of a shorter time constant may attenuate waveforms.
Filter Settings	Low Hz = 0.1 Hz; High Hz = 5 Hz; Gain: Output range is +/- 100 Microvolts.
Sensitivity	Response is dependent on a variety of variables such as Sensor application, patient effort, Probe placement, patient airflow, ambient room temperature, CPAP pressure, etc.

Cannula:

Description	Properties
Cannula	Soft, clear and light for a comfortable patient fit.
Patient attachment	Over-the-ear, "Bolo-Style" adjustable slide.
Length	213 cm (7') kink resistant tubing.
Terminator	Male luer lock.

Caution Statement:

US Federal Law restricts this device to sale by or on the order of a physician. This product is for diagnostic purposes only and is not intended to be used as an apnea monitor or in a life supporting or sustaining capacity.

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©2021 Cadwell Industries, Inc. All rights reserved. PN 302993-639-2 Rev. 00
Cadwell Industries Inc.
1 (800) 245-3001 U.S.A. toll-free | +1 (509) 735-6481
This document contains trademarks that belong to Cadwell Industries, Inc.