
Cadwell® Sleepmate

Instructions for Use (IFU)

Beltless Piezo Respiratory Effort Sensor

Contents**Page**

English	4
Dansk	6
Deutsch	8
Español	10
Suomi	12
Français	14
Italiano	16
Nederlands	18
Norsk	20
Português	22
Svenska	24

Intended use

The Cadwell Sleepmate™ Respiratory Effort Sensor is intended to detect respiratory effort from chest or abdomen movement during sleep research and clinical studies of sleep. This Belt is intended for use in supervised medical settings.

Instruction for use

Place the sensor on the patient's chest or abdomen. Use a piece of adhesive tape over each "wing" of the sensor to attach it to the patient's chest to acquire a strong and stable Effort Signal. It is ok to allow the tape to lap over on to the main sensor body; however do not use a single piece of tape across the entire sensor or allow the 2 pieces of tape to overlap each other. This will greatly reduce the output signal or result in no signal. Secure wire in a manner to avoid entanglement with patient. NO interface or external power source is necessary because the sensor produces its own signal. Connect the sensor to the appropriate channel input of your recorder.

Ensure the sensor is securely attached to the patient as described and you have clear, reliable signals prior to lights out.

Cleaning

You can clean the system by rubbing with any cleansing agent that is not corrosive toward plastic.

The sensor must be completely dry before using again to insure patient safety. For sterilization implement a standard procedure for gas sterilization. Do not autoclave this product.

Use the molded ends of the connectors when removing the sensor from the junction box.

Technical Specifications

Description	Properties
Lead wire	Two conductors white, approximately 80 cm (31.5").
Connectors	Two 1.5 mm (0.1") Safety pins, some models may use special connectors
Signal output	Output signal is typically 1 mV for normal respiration low Frequency
Filter/Time constant	10 Hz / One second or longer. The use of a shorter time constant may attenuate waveforms
High Frequency Filter	0.1 Hz
Output Impedance	75 k Ω + 25 k Ω

Caution Statement:

US Federal Law restricts this device to sale by or on the order of a physician. This product is for diagnostic purposes only and is not intended to be used as an apnea monitor or in a life supporting or sustaining capacity

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©2021 Cadwell Industries, Inc. All rights reserved. PN 302995-639-1 Rev. 00
Cadwell Industries Inc.
1 (800) 245-3001 U.S.A. toll-free | +1 (509) 735-6481
This document contains trademarks that belong to Cadwell Industries, Inc.