
Cadwell® Sleepmate

Instructions for Use (IFU)

Upper Airway Resistance System

Contents**Page**

English	4
Dansk	6
Deutsch	8
Español	11
Suomi	14
Français	16
Italiano	18
Nederlands	21
Norsk	23
Português	25
Svenska	27

Intended use

The Cadwell Sleepmate™ Upper Airway Resistance Sensor System is designed to measure the movement in the sternal notch associated with respiratory effort. The sensor take advantage of a new patent pending technology that produces a small voltage when movement occurs. The system produces its own signal with no danger to the patient. This voltage, when amplified by the physiograph, produces a clean and reliable measurement of upper airway resistance.

Instruction for use

1. **Attach Sensor to Cable:** Connect the sensor to the cable by grasping the sensor at the exposed silver end, and insert at least one half of the sensor's silver end into the interface cable. *Note: even after complete insertion a portion of the connector may remain exposed.*
2. **Attach Sensor to Patient:** Identify the patient's suprasternal notch. Clean the patient's skin around the notch with a gentle nonabrasive cleanser. Peel off the sensor paper backing to expose the adhesive portion. Place the sensor in the suprasternal notch of the patient. The outer rim of the sensor needs to be anchored to the bones that form the frame for the sensor. The signal is achieved by measuring the movement of the sensor, which is adhered to the skin. The greater the upper airway resistance the greater the movement of the sensor. Therefore, it is critical that the sensor be securely attached to the skin via the adhesive disk.
3. **Connect the Cable to a Recorder:** Connect the two lead wire connectors to an active and reference pin in the recorder and select an AC input channel on the recorder. NO interface or external power source is necessary because the sensors produce their own signal.
4. **Ensure the sensor:** Ensure the sensor is securely attached to the patient as described and you have clear, reliable signals prior to lights out.
5. **Removal:** Carefully disengage the sensor from the patient's neck. The sensor was designed for single patient use.

The cable connectors provide a direct connection to the recorder. The sensor fits most patients. However, the soft foam sensor collar can be cut, if necessary, to adapt the sensor to fit all patients.

Cleaning

Please clean the cable by rubbing with any cleansing agent that is not corrosive toward plastic. **The entire assembly must be completely dry before using again to insure patient safety.** For sterilization implement a standard procedure for gas sterilization. Do not autoclave this product. Use the molded ends of the connectors when removing the connector from the junction box.

Technical Specifications

Description	Properties
Cable	Two connectors white lead, approx. 200 cm. (80")
Battery life	Average life with daily use is approx. 1 and 1/2 years
Connectors	1.5 mm (0.1") safety pins.
Signal output	Output signal is typically 0.5 mV PTP, set gain accordingly.
AMP settings	Low frequency = 0.1 Hz; High frequency = 5 Hz.
Sensor	Oval shaped foam disk with laminated adhesive back. Surface area = 60 mm width X 30 mm height (2.4" x 1.2").
Battery	3V Lithium (integrated in the device)

Caution Statement:

US Federal Law restricts this device to sale by or on the order of a physician. This product is for diagnostic purposes only and is not intended to be used as an apnea monitor or in a life supporting or sustaining capacity.

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©2021 Cadwell Industries, Inc. All rights reserved. PN 302995-639-5 Rev. 00
Cadwell Industries Inc.
1 (800) 245-3001 U.S.A. toll-free | +1 (509) 735-6481
This document contains trademarks that belong to Cadwell Industries, Inc.